

Perfect Bleach: Effective and gentle



Simple and user friendly

The application of Perfect Bleach is uncomplicated and user friendly: the patient simply fills the custom-fabricated and optimal fitting trays with the high viscosity gel and wears them as instructed. Wearing the trays for two hours is sufficient for the gel to achieve full effect.

Initial results are already evident after the second application. As the

viscosity of the gel ensures it does not run out of the tray, gingival irritation is prevented. The fresh mint flavour is also pleasant for the patient. Furthermore, the concentration of the carbamide peroxide gel remains stable even if the cold chain is interrupted.

Complete Set and Patient Kit

The Perfect Bleach Complete Set contains all the components needed: six 2.4 ml syringes of whitening gel, two deep-drawing films, one 1.2 ml syringe of Block Out Gel LC for the fabrication of the trays, application pictogram card, shade guide to verify treatment success, storage case for the trays, and a cosmetic bag. The Patient Kit contains four 2.4 ml syringes of whitening gel, application pictogram card, shade guide, storage case for the trays, and a cosmetic bag.

With Perfect Bleach, VOCO offers a home whitening system that meets all of the requirements of modern tooth whitening and thus the highest aesthetic demands in an equally simple and cost-effective way. Vital teeth that are discoloured because of ageing or diet-related accretion can be whitened quickly and effectively. Devitalised teeth can also be whitened externally and internally. The treatment of tetracycline discolouration, superficial enamel discolouration from fluorosis, and discolouration caused by trauma due to bleeding is also possible. Perfect Bleach is a carbamide peroxide-based home whitening material. Depending on the type of discolouration, either a 10 % or 17 % concentration can be used.

Gentle and effective

Distinct whitening success can be achieved with whitening gels containing low concentrations of carbamide peroxide and with a simultaneous gentle treatment of the tooth substance. The carbamide peroxide concentration in Perfect Bleach is sufficiently high to provide effective and gentle treatment of discolouration. The sodium fluoride and potassium in 17 % Perfect Bleach ensure that hypersensitivity is prevented during and after the treatment. Perfect Bleach provides a thorough oxidation of accumulated pigments without damaging the tooth substance or existing restorations.

contact

cosmetic
dentistry

VOCO GmbH
P.O. Box 767
27457 Cuxhaven
Germany

Website: www.voco.com