


Straumann® Pro Arch BLX

New confidence, flexibility and efficiency with one system.

One system for most full-arch challenges of the immediate rehabilitation with the Straumann® BLX Implant System designed for primary stability.



CONFIDENCE

Implant design and Dynamic Bone Management™ for primary stability, Roxolid® for strength^{1,2} and SLActive® for faster healing³.



FLEXIBILITY

Broad range of 42 implants and optimized shallow threads of narrow implants designed for use in the soft and hard bone.



EFFICIENCY

One connection for all diameters. Ø 3.75 mm implant for all indications and slim and smooth Screw-retained Abutments.

REFERENCES

1 Steinemann S.G. 'Titanium - the materials of choice?' Periodontology 2000, Vol. 17, 1998, 7-21 2 Data on file for Straumann cold-worked titanium and Roxolid® Implants SLActive 3 Lang NP, Salvi GE, Huynh-Ba G, Ivanovski S, Donos N, Bosshardt DD. Early osseointegration to hydrophilic and hydrophobic implant surfaces in humans. Clin Oral Implants Res. 2011 A



Experience **Straumann® Pro Arch BLX**.
Contact your local Straumann
representative or scan QR-Code



Supported by digital service **Smile in a Box™**

A0026/en/B/00-IIM 10/22