

initial™ Zirconia Coloring Liquid

Customised, your way!

Faster, more predictable
and highly individualised
following the Initial IQ philosophy

Smaller inventory, greater gains
Obtain every shade fast and easy
with a reduced inventory



Brush or dip
Two ways, your call

The **Initial Zirconia Coloring Liquid** system consists of 12 water-based colouring liquids and 4 marking colours following the 'IQ philosophy' for the replication of natural tooth colours and effects prior to sintering. Simply apply the liquids onto the monolithic or anatomically reduced restoration after milling using a painting or dipping technique.

