

# CLEARFIL S<sup>3</sup> BOND PLUS

## High-tech and high speed for high performance



**CLEARFIL S<sup>3</sup> BOND PLUS** is the one-step adhesive for direct restorations. It has been developed based on CLEARFIL S<sup>3</sup> BOND technology, which reveals five years of clinical evidence. Etching, priming and bonding steps are completed with one-liquid, one-coat treatment.

CLEARFIL S<sup>3</sup> BOND PLUS is easy and comfortable to use because time-intensive work steps, such as shaking the bottle, exact mixing and application of several components, multiple layering and rubbing in on the tooth surface, are no longer necessary.

The high bond strength of CLEARFIL S<sup>3</sup> BOND PLUS is the result of two Kuraray innovative technologies: the high performance initiator and the proven MDP monomer. The new high performance initiator for light-curing improves clinical performance by building up more active radicals than conventional initiators and making the bond more impervious to water thanks to the increased polymerisation ratio.

The excellent bond strength to enamel and dentine of CLEARFIL S<sup>3</sup> BOND PLUS marks a new standard among one-step adhesives, enabled by Kuraray's innovative adhesive technology, which has been successful for more than 30 years. Its fast and simple applica-

The integrated adhesive monomer MDP ensures a strong chemical bond to hydroxyapatite. In use for more than 25 years in successful products like PANAVIA F 2.0 and CLEARFIL SE BOND, the MDP monomer attains a high bond strength and reliable adhesive durability to tooth structure.

tion combined with its fluoride-releasing property makes CLEARFIL S<sup>3</sup> BOND PLUS the ideal choice for all kinds of clinical situations, including paediatric therapy and core build-ups with the new CLEARFIL DC CORE PLUS. Furthermore, the system allows for precise and economic dosing.

The stable interface between tooth and composite resin is furthermore ensured by Kuraray's unique Molecular Dispersion Technology contained in CLEARFIL S<sup>3</sup> BOND PLUS. The technology allows hydrophilic and hydrophobic components of the bond to be homogeneously combined—constantly. Thus, there is no phase separation, resulting in the formation of a homogenous layer on the preparation surface. The bond, showing no water voids, excels through its excellent bond performance.

CLEARFIL S<sup>3</sup> BOND PLUS is always ready to use and will have users persuaded of its durable high bond strength.

The main benefit of using one-step adhesives is to reduce technical errors, thanks to a short application time and simple handling. The time-saving procedure makes CLEARFIL S<sup>3</sup> BOND PLUS ready for immediate use in only three short steps: apply, air-dry, light-cure—all done in less than 30 seconds.



_contact	<b>cosmetic</b> dentistry
<p><b>Kuraray Europe GmbH</b>                  Industriepark Höchst, Building F 821                  65926 Frankfurt/Main                  Germany</p> <p>dental@kuraray.eu                  www.kuraray-dental.eu (EU)                  www.kuraraydental.com (USA)</p>	



**APRIL 26- 28, 2012 • BELLA CENTER • COPENHAGEN • DENMARK**

**Welcome to the 45th Scandinavian Dental Fair**  
**The leading annual dental fair in Scandinavia**

# SCANDEFA 2012



**The 45th SCANDEFA invites you to exquisitely meet the Scandinavian dental market and sales partners from all over the world in springtime in wonderful Copenhagen**

## **SCANDEFA 2012**

Is organized by Bella Center and is being held in conjunction with the Annual Scientific Meeting, organized by the Danish Dental Association ([www.tandlaegeforeningen.dk](http://www.tandlaegeforeningen.dk)).

175 exhibitors and 11.422 visitors participated at SCANDEFA 2011 on 14,220 m<sup>2</sup> of exhibition space.

## **Exhibit at Scandefa**

Book online at [www.scandefa.dk](http://www.scandefa.dk)  
Sales and Project Manager, Christian Olrik  
[col@bellacenter.dk](mailto:col@bellacenter.dk), T +45 32 47 21 25

## **Travel information**

Bella Center is located just a 10 minute taxi drive from Copenhagen Airport. A regional train runs from the airport to Orestad Station, only 15 minutes drive.

## **Check in at Bella Center's newly built hotel**

Bella Sky Comwell is Scandinavia's largest design hotel. The hotel is an integral part of Bella Center and has direct access to Scandefa. Book your stay on [www.bellasky.dk](http://www.bellasky.dk)