

RECIPROC wins innovation prize 2011



_RECIPROC is probably the most significant new development since nickel-titanium instruments were first introduced to rotary preparation of root canals. RECIPROC completely prepares and shapes the root canal with one single instrument. The system was developed to simplify the procedure of preparing root canals while ensuring maximum security during the process. Hence, since its launch in March 2011, the RECIPROC® system has been adopted in numerous dental offices and clinics replacing manual instrumentation as well as rotary NiTi systems.

At six German dental trade shows, dentists voted for RECIPROC one-file endo to be awarded first prize in the materials and instruments for dental treatment

category of DZW and pluradent's German INNOVATIONSPREIS 2011.

During the ADF congress in Paris (November 2011), out of 56 entries Réciprocité (reciprocation technique) won the Dental Profession's Special Innovation Award.

It is the simplicity of working without instrument changes and only one motor setting, combined with safety and economic single use, that has dentists convinced.

_contact	roots
VDW Fax: +49 89 62734-304 info@vdw-dental.com www.reciproc.com	

Everything you like about X-Smart— with a Plus

_The new generation of the popular X-Smart endo motor from DENTSPLY Maillefer, the X-Smart Plus, will allow you to discover the WaveOne reciprocating single-file technique advocated by leading endodontists. The X-Smart user interface has been further improved by a large, bright colour screen, with a colour-coded file library for file selection at a single glance, making it the endo motor of choice for all ProTaper Universal users. The X-Smart Plus retains the familiar X-Smart features such as the miniature contra-angle head and the on/off button on the motor hand-piece.

as well as eight free programmes for individual settings. It is operated by a rechargeable battery and in continuous rotation provides a speed range of between 250 and 1,200 rpm and a torque range of between 0.6 and 4 Ncm. A warning sound helps you keep track of the file stress and the auto-reverse rotation at the torque limit reduces the risk of file breakage.

The X-Smart Plus endo motor from DENTSPLY Maillefer comes with a three-year warranty.



_contact	roots
DENTSPLY Maillefer www.dentsplymaillefer.com	

The motor works in both reciprocating motion and continuous rotation. The file library contains preprogrammed settings for WaveOne, ProTaper Universal, PathFile, Gates and RECIPROC,

implants

international magazine of oral implantology



You can also subscribe via
www.oemus.com/abo



Subscribe now!

One issue free of charge!

☐ I hereby agree to receive a free trial subscription of **implants international magazine of oral implantology** (4 issues per year).

I would like to subscribe to **implants** for € 44 including shipping and VAT for German customers, € 46 including shipping and VAT for customers outside of Germany, unless a written cancellation is sent within 14 days of the receipt of the trial subscription. The subscription will be renewed automatically every year until a written cancellation is sent to OEMUS MEDIA AG, Holbeinstr. 29, 04229 Leipzig, Germany, six weeks prior to the renewal date.

▶ **Reply** via Fax +49 341 48474-290 to OEMUS MEDIA AG or via E-mail to grasse@oemus-media.de

Last Name, First Name

Company

Street

ZIP/City/County

E-mail

Signature

Notice of revocation: I am able to revoke the subscription within 14 days after my order by sending a written cancellation to OEMUS MEDIA AG, Holbeinstr. 29, 04229 Leipzig, Germany.

Signature

OEMUS MEDIA AG
Holbeinstr. 29, 04229 Leipzig, Germany, Tel.: +49 341 48474-0, Fax: +49 341 48474-290, E-Mail: grasse@oemus-media.de



Moisture-activated temporary filling and sealing material

Tempit is packaged in 30 tips of 0.35g each in a jar.

_contact

roots

Centrix

Tel.: +49 221 530978-0

koeln@centrixdental.com
www.centrixdental.com

_Still squeezing those messy tubes? With Tempit, our unit-dose, moisture-activated temporary filling and sealing material, simply inject into a moist prep, tamp down and then dismiss the patient. Tempit sets in five minutes or less.

Tempit was developed specifically for sealing access cavities between root-canal treatments. It seals medicaments in and contaminants out between root-canal treatments.



Pre-filled, it eliminates set-up time. Unique delivery provides:

- _direct injection without worrying that the material will harden before use;
- _no clean-up; and
- _eliminates the voids you get when placing other materials.

No more spatulas! No more tubes or jars!_

AD

cosmetic dentistry _ beauty & science

One issue free of charge!



Subscribe now!

You can also subscribe via
www.oemus.com/abo

☐ I hereby agree to receive a free trial subscription of **cosmetic dentistry** (4 issues per year). I would like to subscribe to **cosmetic dentistry** for € 44 including shipping and VAT for German customers, € 46 including shipping and VAT for customers outside of Germany, unless a written cancellation is sent within 14 days of the receipt of the trial subscription. The subscription will be renewed automatically every year until a written cancellation is sent to OEMUS MEDIA AG, Holbeinstr. 29, 04229 Leipzig, Germany, six weeks prior to the renewal date.

► **Reply** via Fax +49 341 48474-290 to OEMUS MEDIA AG or per E-mail to grasse@oemus-media.de

Last Name, First Name

Company

Street

ZIP/City/Country

E-mail

Signature

Notice of revocation: I am able to revoke the subscription within 14 days after my order by sending a written cancellation to OEMUS MEDIA AG, Holbeinstr. 29, 04229 Leipzig, Germany.

Signature

roots 2/12





RECIPROC wins innovation prize 2011



RECIPROC is probably the most significant new development since nickel-titanium instruments were first introduced to rotary preparation of root canals. RECIPROC completely prepares and shapes the root canal with one single instrument. The system was developed to simplify the procedure of preparing root canals while ensuring maximum security during the process. Hence, since its launch in March 2011, the RECIPROC® system has been adopted in numerous dental offices and clinics replacing manual instrumentation as well as rotary NiTi systems.

At six German dental trade shows, dentists voted for RECIPROC one-file endo to be awarded first prize in the materials and instruments for dental treatment category of DZW and pluradent's German INNOVATIONS PREIS 2011.

During the ADF congress in Paris (November 2011), out of 56 entries Réciprocité (reciprocation tech-

nique) won the Dental Profession's Special Innovation Award.

It is the simplicity of working without instrument changes and only one motor setting, combined with safety and economic single use, that has dentists convinced.

_contact	roots
VDW Fax: +49 89 6273 4304 info@vdw-dental.com www.reciproc.com	

Dental Tribune pour iPad – Votre sélection de nouvelles hebdomadaires

Chaque semaine, nos éditeurs choisissent les meilleurs articles et vidéos des différentes éditions des Dental Tribune dans le monde.

Créez votre édition personnelle, dans la langue de votre choix, en sélectionnant uniquement les publications et vidéos qui vous intéressent.

ipad.dental-tribune.com

