

Moisture-activated temporary filling and sealing material

Tempit is packaged in 30 tips of 0.35g each in a jar.

Still squeezing those messy tubes? With Tempit, our unit-dose, moisture-activated temporary filling and sealing material, simply inject into a moist prep, tamp down and then dismiss the patient. Tempit sets in five minutes or less.

Tempit was developed specifically for sealing access cavities between root-canal treatments. It seals medicaments in and contaminants out between root-canal treatments.



Pre-filled, it eliminates set-up time. Unique delivery provides:

- direct injection without worrying that the material will harden before use;
- no clean-up; and
- eliminates the voids you get when placing other materials.

No more spatulas! No more tubes or jars!

_contact roots

Centrix

Tel.: +49 221 530978-0

koeln@centrixdental.com
www.centrixdental.com

RAYPEX 6— State-of-the-art technology

For the precise location of the apical constriction, electronic means have been proven to be superior to using radiographs. RAYPEX devices ensure precise measurement and feature a unique zoom-in function, which significantly aids endodontic treatment. The latest model, RAYPEX 6, redefines user-friendliness and design standards.

Its design, operation and user-friendliness were developed to a high standard, resulting in a stylish and handy casing, which can be folded flat when not in use. The high-resolution touch-screen with a self-explanatory user-interface is comparable with a smartphone, and the menu guidance was designed to ensure precise and fast navigation through the functions.

For precise length determination, RAYPEX 6 uses the latest multi-frequency technology. The 3-D graphic display constantly indicates the

position of the tip of the file, while the automatic zoom-in function in the region of the apical constriction enhances safety in the critical section just before reaching the working length. Acoustic signals allow for a "blind" localisation of the apex.

RAYPEX 6 is ready for use—simply plug-in and start. Any personalisation of display colour and tones, for instance, can be done quickly.

RAYPEX 6 means precise length determination and touch-screen convenience.



_contact roots


VDW GmbH

Bayerwaldstr. 15
81737 Munich
Germany

info@vdw-dental.com
www.vdw-dental.com

implants

international magazine of oral implantology



You can also subscribe via
www.oemus.com/abo

 **Subscribe now!**

I hereby agree to receive a free trial subscription of **implants international magazine of oral implantology** (4 issues per year).

I would like to subscribe to **implants** for € 44 including shipping and VAT for German customers, € 46 including shipping and VAT for customers outside of Germany, unless a written cancellation is sent within 14 days of the receipt of the trial subscription. The subscription will be renewed automatically every year until a written cancellation is sent to OEMUS MEDIA AG, Holbeinstr. 29, 04229 Leipzig, Germany, six weeks prior to the renewal date.

▶ **Reply** via Fax +49 341 48474-290 to OEMUS MEDIA AG or via E-mail to grasse@oemus-media.de

One issue free of charge!

Last Name, First Name	
Company	
Street	
ZIP/City/County	
E-mail	Signature

Notice of revocation: I am able to revoke the subscription within 14 days after my order by sending a written cancellation to OEMUS MEDIA AG, Holbeinstr. 29, 04229 Leipzig, Germany.

Signature

