


Want to save some space?

Cut the excess with our **universal adhesive solutions**


Restore

- ✓ Compatible with all GC composites
- ✓ High bond strength in all etching modes




Lute

- ✓ Compatible with all GC resin cements
- ✓ Very low film thickness




Build up cores


- ✓ Optional dual cure mode available
- ✓ High bond strength, even in bulk-fill application



Repair

- ✓ Suitable for direct and indirect restorations
- ✓ Repair in simple steps





G-Premio BOND

- ✓ Clean application
- ✓ 300 drops/bottle
- ✓ Universal applicability



G-Premio BOND

One component light-cured universal adhesive



GC composite range

light-cured restorative composites



G-CEM ONE

Universal dual-cured composite resin cement



GRADIA CORE

Dual-cure composite for core build-up & post cementation



everX Flow

Short-fibre reinforced flowable composite



Repair kit

Chairside kit for all intraoral repairs

Solved with GC!

