

# Digitally networked implant treatment

NSK has developed and improved the Surgic Pro surgical micro-motor system, which is now available as Surgic Pro2. Surgic Pro2 provides a Bluetooth connection to the Osseo 100+ osseo-integration gauge and a Link Set connection to the VarioSurg3 ultrasonic surgical system. By installing an application and connecting an iPad to the control unit, detailed intervention data can be displayed in real time. Both Surgic Pro2 treatment parameters and Osseo 100+ data can be displayed and stored on the iPad. The traceability of patient-specific treatment data means that implant treatment can be customised for each patient. The Surgic Pro2 and VarioSurg3 can be operated wirelessly and hygienically using a common foot control. But even without connection to other units, Surgic Pro2 offers advantages for the practitioner. The micromotor is short, lightweight and ergonomic. A high-resolution LED in the micromotor serves as light source to match working conditions under full daylight. The pump module is quiet, and the irrigation tube is easy to attach. The Surgic Pro2 is simple and intuitive to use, making it easy to upgrade from the previous model to Surgic Pro2.

Take advantage of the current savings packages in the Get it! campaign:



## Contact address

NSK Europe GmbH  
+49 6196 77606-0  
[www.nsk-europe.de](http://www.nsk-europe.de)

\* The articles in this category are provided by the manufacturers or distributors and do not reflect the opinion of the editorial team.