

Success made easy for everyone



91%
clinical
success
after 2 years
—
Exceeds
success rates found
in the literature*

Your Bioactive Root Canal Sealer



- Ready-to-use hydraulic sealer in a syringe
- No shrinkage - bioactive - biocompatible - high pH
- Whatever your technique, cold or warm

*Clinicaltrials.gov: Root canal obturation with a ready-to-use root canal sealer (PA1704) versus BioRoot™ RCS: a randomized controlled trial.

Class III Medical Device - Certified by BSI (2797) for MDR/EU compliance.

Please visit our website for more information
www.septodont.com

Visit us
HALL 4.1

