

An implant like no other.

High initial stability even in compromised bone situations.

Maximum bone preservation.

Excellent esthetics.

PREVIEW

Ø 3.0 mm for limited spaces
Launch in fall 2011



It's called NobelActive. The drilling procedure of this unique implant is designed to retain as much bone as possible, while the implant body and thread design condenses bone during insertion enhancing initial stability. The back-tapered coronal design and built-in platform shifting maximize alveolar bone and soft tissue volume, and the sharp apex and cutting blades enable you to adjust

the implant position for optimal restorative orientation. Together with the strong conical connection and a comprehensive prosthetic assortment, NobelActive allows you to achieve excellent esthetic results. After 45 years as a dental innovator we have the experience to bring you future-proof and reliable technologies for effective patient treatment. **Their smile, your skill, our solutions.**



Visit nobelbiocare.com/nobelactive

