

**Precision demands
the right equipment.**

CLEARFIL™ SE BOND

The adhesive that meets your high demands – with a durable marginal sealing thanks to the exceptionally high hydrolysis resistance & with least postoperative sensitivity.

CLEARFIL™ SE BOND is Kuraray's high performance, self-etching adhesive for direct restorations – known for its outstanding long-term marginal sealing thanks to the product's powerful hydrolysis resistance. Furthermore, dentists and patients value the product's least postoperative sensitivity.

CLEARFIL™ SE BOND enables an excellent long-term performance due to the high sealing ability of the bonding agent and the mild-etching, low technique sensitive self-etching primer.



With the in vitro study on the marginal adaptation to enamel and dentin*¹, CLEARFIL™ SE BOND shows its effectiveness by an outstanding, durable marginal adaptation. Moreover, in 2010, CLEARFIL™ SE BOND proves its exceptional strength within an 8 year clinical study*² showing a retention rate of 97%.

Thus, CLEARFIL™ SE BOND is appreciated as the Gold Standard amongst several universities and opinion leaders.

*¹ Source: K. Huber, G. C. Lier, B. Bott, and M. Hanning: Marginal Adaptation of Composite Resin Restorations Using Self-etching Adhesives, Joint Meeting of the Continental European, Israeli, and Scandinavian (NOF) Divisions of the IADR, Abstract #59, August 25-28, 2004.

*² Source: B. Van Meerbeek, M. Peumans, A. Poitevin, A. Mine, A. Van Ende, A. Neves, and J. De Munck: Relationship between bond-strength tests and clinical outcomes, Dental Materials 26, e100–e121, 2010.