Vol. 6 • Issue 4/2012

dentistr uty & science

2012

I CE article
The digitized occlusion:

Using something old with something new

Cost effectiveness in implant dentistry

The same-day tooth:

From the diagnosis to the final restoration

०२७११।उ

kuraray



SA CEMENT Automix by Kuraray has now been relaunched as part of the PANAVIA $^{\text{TM}}$ brand with an expanded range of indications including cementation of restorations on implant abutments.

PANAVIA™ SA CEMENT Automix is a dual-curing, self-adhesive resin cement in an auto-mix syringe for routine restorations - the guarantee for you and your patient's satisfaction. The easy removal of the excess paste reduces the risk of injuring the gingival sulcus and damages of the prosthetic.

This is due to Kuraray's adhesive technology based on a proprietary adhesive phosphate monomer "MDP". This unique adhesive monomer from Kuraray provides high bond strength and secures low technique sensivity for consistent, reliable results. The high mechanical stability gives excellent marginal seal for long lasting restorations.

The PANAVIA™ brand can look back on a scientific and clinical track record spanning more than 25 years!