



page 06



page 10



page 26

| editorial

03 **Cologne** in the eyes of **laser dentistry**
| Dr Georg Bach

| overview

06 **Laser versus** conventional therapies
| Cristiane Meira Assunção *et al.*

| research

10 **Photodynamic** therapies – **Blue** versus **Green**
| Dr Michael Hopp *et al.*

| case report

26 **Becoming kissable:**
Laser-assisted haemangioma removal
| Dr Darius Moghtader

| industry report

30 **Er:YAG Garnet** in laser-assisted **crown lengthening**
| Dr Avi Reyhanian

34 **Long-term treatment of peri-implant lesions in geriatric dentistry**
| Dr Georg Bach

38 **X-Runner Er:YAG dental laser** application
| Prof. Dr Carlo Fornaini

| education

44 **1st International Congress of WALED and GLOBAL 2013 in Istanbul**
| AALZ Germany

| meetings

46 **The next chapter in the IDS success story**
| Koelnmesse

| news

42 **Manufacturer News**
48 **News**

| about the publisher

50 | imprint



Cover image courtesy of Fotona, www.fotona.com



page 34



page 44



page 46



After endodontic laser treatment there is no smear layer around the opening of the lateral canal.

The universe at your fingertips.



LightWalker®

Introducing the highest technology dental laser system

Supreme clinical results in:

- TwinLight™ Perio Treatments (TPT)
- TwinLight™ Endo Treatments (TET)
- Effective oral implant surgery
- Gentle TouchWhite™ bleaching
- Patient-friendly conservative dentistry

Unmatched simplicity of use:

- Pre-sets for over 40 applications
- Intuitive user navigation
- Balanced and weightless OPTOflex arm
- Nd:YAG handpiece detection system
- Er:YAG scanner ready

Journey into a new dental experience with speed, precision and great results.
Visit www.lightwalkerlaser.com today!

