




| editorial

- 03 Does **biology still** matter?
Prof. (CAI) Dr Roland Hille

| research

- 06 **Immediate** restoration in the **digital** workflow
José Eduardo Maté Sánchez de Val &
José Luis Calvo Guirado

| case report

- 12 Natural **preservation** of the **emergence** profile
Dr Umut Baysal & Dr Arzu Tuna
- 18 Immediate **temporisation** using **tapered** implants
Dr Mohamad El Moheb

| industry

- 22 Overcoming a **difficult aesthetic** situation
Prof. Dr Dr M. Kemal Unsal
- 28 From **granules** to **foam**
Prof. Dr Dr Stefan Schermer

- 34 **Xenogeneic** bone **grafting** materials
Dr Mike Barbeck, Dr Ronald Unger,
Prof. Dr Frank Witte, Prof. Dr Sabine Wenisch &
Prof. Dr Dr Reiner Schnettler

- 38 **Dynamic** navigation in **fully** edentulous **maxilla**
Prof. Dr Hakan Uysal & Dr Noyan Başal

| events

- 52 Arabian **flair** and **French** charm
Dr Rolf Vollmer
- 54 Joint **meeting** of the **EAO** and **SEPES**

| news

- 42 manufacturer **news**
- 56 **news**

| about the publisher

- 58 imprint



Cover image courtesy of bredent medical GmbH & Co. KG
www.bredent-medical.com

a perfect fit™

camlog

THE NEW COMFORT CLASS.

The COMFOUR® Abutment System for occlusal screw-retained restorations is an ideal addition to the CAMLOG® and the CONELOG® Implant System. With the COMFOUR® System, you can offer your edentulous patients the option of an immediate, comfortable and permanent denture on four or six implants – and thus a considerable gain in quality of life. www.camlog.com

CAMLOG®
SYSTEM

CONELOG®
SYSTEM



Aligning tool



Scanning cap



Titanium cap



Flexible handle

COMFOUR® offers you a range of options to overcome the challenges you face in your implant practice, making them easier while saving time.